

Technical Evaluation of **HPGP & HDROZONE** (Hydronium Heavy Molecule) Sterilization Methods in Flexible Endoscopes

Hasan Tahsin ÖZBEK • Bülent DEVECİ • Yüksel ERGÜN

White Paper - **24.09.2025**

SUBJECT

This document evaluates the technical suitability, risk areas, and overall outcomes of VH_2O_2 , HPGP, and HDROZONE sterilization approaches for flexible endoscopes, taking into account the material structure of the devices.

The evaluation is structured around the following key aspects:


- Endoscope material structure
- Critical surfaces in sterilization
- Intra-lumen sterilization requirements
- Behavior of polymer, elastomer, metal, and optical components
- Suitability and risk assessment for HDROZONE




KEYS

HPGP, VH_2O_2 , Hydronium, Infiltration, Diffusion, Condensation, Hybrid, Distal

1. WHAT IS THE PROBLEM IN ENDOSCOPE STERILIZATION?




ENDOSCOPE CROSS SECTION



Flexible endoscope sterilization faces a fundamental challenge: the device is not made of a single homogeneous material.

A flexible endoscope is a complex multi-layered hybrid medical device that includes:

- Outer sheath
- Inner spiral / braid structure
- Deflectable distal tip
- Working channel / lumens
- Optical system
- Control cables
- Adhesives and connectors



Because of this structure, sterilization is not merely an external surface issue.

1.1 Approach of H_2O_2 (VH_2O_2) and HPGP Sterilization Methods to These Problems


- Vapor-phase hydrogen peroxide creates strong oxidative effect on surfaces.
- It relies on penetration, diffusion, and condensation mechanisms.
- Reactive loss can occur inside lumens.
- Efficacy against biofilm depends on surface access.

CONCLUSION:

H_2O_2 systems are strong in surface sterilization but show limited effectiveness in long and narrow lumens.

2. MATERIAL-BASED PROBLEM DEFINITION IN ENDOSCOPES

2.1 Outer Sheath



Materials used:

- TPU
- Silicone elastomer
- Pebax



Conclusion:

The outer sheath is the primary critical area for sterilization-related damage.

H₂O₂ vs Hdrozone

- H₂O₂:
 - High oxidative density, increased risk of surface damage
- Hdrozone:
 - Controlled oxidant distribution, low-temperature advantage


2.2 Inner Spiral / Braid Structure



CONCLUSION:

The metal carrier structure is a low-risk area for both HDROZONE and H₂O₂.

2.3 Distal Deflectable Tip





CONCLUSION:

The distal tip is a high-risk area where chemical and mechanical stresses converge.


H₂O₂ vs Hdrozone

- H₂O₂: Possible condensation and localized concentration increase
- Hdrozone: More homogeneous gas-phase distribution


2.6 Adhesives and Bonding Agents



a) Adhesive or interface failure



b) Cohesive failure



c) Substrate or adherend failure




CONCLUSION:

Epoxy and adhesives are among the most critical weak points in endoscopes.

H₂O₂ vs Hdrozone

- H₂O₂: Intense oxidative effect
- Hdrozone: Controlled chemical effect

TECHNICAL SUITABILITY ASSESSMENT OF H₂O₂ VS. HDROZONE IN ENDOSCOPE STERILIZATION



REGION	H ₂ O ₂	HDROZONE
Outer Sheath	Medium-High Risk	Medium Risk
Metal Structure	Low Risk	Low Risk
Distal Tip	High Risk	Medium-High Risk
Lumen	Limited	Higher Potential
Epoxy	High Risk	Medium Risk

Key Determining Factors for HDROZONE:

- Lumen inner surfaces
- Distal tip
- Biofilm regions


Comparative Result

- H₂O₂: Surface-dominant sterilization
- Hdrozone: Lumen-focused sterilization approach

MATERIAL DAMAGE COMPARISON (H₂O₂ VS. HDROZONE)

- H₂O₂: High oxidative density
- Hdrozone: More balanced oxidant distribution

EFFICACY COMPARISON OF H₂O₂ VS. HDROZONE IN ENDOSCOPE STERILIZATION




- H₂O₂: Strong surface effect, low penetration potential
- Hdrozone: Strong surface effect, high penetration potential

ULTRA-LOW TEMPERATURE EFFECT IN STERILIZATION

- H₂O₂: Generally 45–55°C (low temperature), but local heating may occur
- Hdrozone: Controlled ultra-low temperature (37–42°C)

LUMEN PERFORMANCE COMPARISON IN ENDOSCOPES



Lumen Sterilization Comparison

- H₂O₂: Limited in long lumens (effective up to ~1.5 m)
- Hdrozone: Gas-phase advantage (effective up to 15 m)

OVERALL COMPARATIVE EVALUATION OF H₂O₂ VS. HDROZONE IN ENDOSCOPE STERILIZATION

Strengths of Both Methods

H₂O₂:

- Strong surface sterilization
- Widespread use and acceptance

Hdrozone:

- Ultra-low temperature
- Superior gas-phase penetration
- Significant lumen advantage

Critical Risks of Both Methods

H₂O₂:

- Epoxy damage
- Excessive surface oxidation

Hdrozone:

- Requires careful process control

FINAL CONCLUSION

- H₂O₂: Surface-oriented sterilization with limited lumen efficacy due to condensation issues.
- Hdrozone: Lumen-oriented hybrid hydronium sterilization method thanks to its hybrid nature.

Overall Assessment:

In flexible endoscope sterilization, H₂O₂ systems are strong in surface sterilization, whereas **HDROZONE systems offer significant technical advantages** in both surface sterilization and long-lumen penetration.

HDROZONE is a hybrid heavy-molecule sterilization method that provides clear superiority over H₂O₂ systems — particularly in lumen penetration and low-temperature compatibility — for flexible endoscopes. It serves as an effective problem-solving technology for complex, long-lumen medical devices.

- 30- Reference The literature sources we used in our hdrOzone studies