

MICROBIOLOGICAL EVALUATION OF THE LOW-TEMPERATURE STERILIZATION METHOD OPERATING WITH HYDRONIUM

Duygu Perçin Renders ¹ , Hasan Tahsin Özbek ²

1.Kütahya Sağlık Bilimleri Üniversitesi Tıp Fakültesi Tıbbi Mikrobiyoloji AD, KÜTAHYA

2. TEKNO MAR Makina İmalat İthalat İhracat Sanayi ve Ticaret Ltd. Şti., ANKARA

INTRODUCTION AND PURPOSE

- Low-temperature sterilization methods are frequently used in Central Sterilization Units (CSUs).

However, they have certain disadvantages, such as:

- Limitations in lumen diameter and length
- Toxicity
- This study was conducted to clinically evaluate,
- The HRF 3000 HdrOZON (Teknomar, Türkiye) device, which is claimed to eliminate these disadvantages,
- Which performs sterilization by using Hydronium ion technology obtained by the combined use of Hydrogen Peroxide and Ozone gases.


MATERIALS AND METHODS



Operating stages of the device

After the chamber temperature reaches 44 °C:


- 1- Vacuum phase:** Bringing the sterilization chamber to -600, -900 mBar vacuum,
- 2- Injection phase:** Vaporization and injection of 60% H₂O₂ concentration,
- 3- Diffusion phase:** Diffusion of hydrogen peroxide into the materials to be sterilized,
- 4- Ozonation:** Taking oxygen into the system, ozonating it and delivering it to the chamber,
- 5- Adsorption:** Adsorption of hydrogen peroxide and ozone onto the materials to be sterilized by forming a heavy molecule.
- 6- Dilution:** Dilution of Hydronium
- 7- Plasma phase:** Performing HRF COLD PLASMA under -300/-700 mBar atmospheric vacuum
- 8- Ventilation:** Aeration phase: return to normal atmospheric pressure.




Sterilization Mechanism

Activation mechanism between H_2O_2 and ozone


Reaction of Ozone with Hydrogen Peroxide (Perozone Process): A Revision of Current Mechanistic Concepts Based on Thermokinetic and Quantum-Chemical Considerations. Activation mechanism of H_2O_2 by ozone, (DORÉ, 1989).



With the application of H_2O_2 and O_3 together, a sterilization effect is achieved on microorganisms through adsorption by the action of hydronium heavy molecules, radicals, hydroxyls, and ions such as HO_3 , HO_4 , H_2O_3 , H_2O_4 , HO_5 , OOH , H_3O , etc.



The mechanism for adsorption on catalyst and oxidation by ozone or HO radical of adsorbed organic



The possible mechanism for HO or other radical species generation by reaction of ozone with reduced metal of catalyst.

Polytetrafluoroethylene (PTFE) and steel lumen tubes were used as the challenging test system (PCD) in the study.



Challenging test systems used


PCD	Length (mm)	Lumen Diameter(mm)
Teflon PTFE	850	1
Teflon PTFE	1200	2
Teflon PTFE	900	0,4
Teflon PTFE	7500	2
Teflon PTFE	10000	2
Teflon PTFE	15000	2
Steel	500	0,7

BIONOVA brand bacterial spore strips containing 10^6 G. *stearothermophilus* were placed at the closed end of the challenging test systems and packaged with 4A brand TYVEK. All test systems were placed in different localizations in each cycle and tested as full and half cycles.




Incubation in Tryptic Soy liquid medium at 55 °C for 7 days.

Full Cycle



Half Cycle



Calculation of the logarithmic reduction factor

$$\text{LRF} = \log_{10} [\text{control}] - \log_{10} [\text{test}]$$

FINDINGS

- In all test systems and in all cycles, the reduction factor was calculated as 8 log compared to the positive control.
- No microorganism growth was observed in any of the Bionova H₂O₂ biological indicators placed inside 23 PCD lumens in each cycle.



CONCLUSION

In all test systems and in all cycles, the reduction factor was calculated as 8 log compared to the positive control. No microorganism growth was observed in any of the Bionova H₂O₂ biological indicators placed inside 23 PCD lumens in each cycle.

THANK YOU